

Zibo Combi Active

Comfort and flexibility - with or without a pump

ZiboCare has developed a new mattress that is the result of 12 years of experience. Introducing the Combi Active Hybrid mattress.

This mattress was developed to address the needs of the broadest possible target group using only the highest quality components and with features focusing on recycling and a better economy for both municipalities and hospitals. For the less mobile and more vulnerable users, a pump can be connected for extra pressure relief using alternating pressure levels.



- **Clinical features:**

Optimal pressure distribution and pressure relief with multiple treatment options. Reduces heat and moisture. Built-in heel function. Transverse air cells provide maximum stability and comfort.

- **The Combi Active is ideal for:**

The treatment of category 1-4 pressure injuries. Prevention for high-risk users, users who cannot / will not reposition themselves, pain users experiencing high levels of pain, sensitive users, people suffering from dementia.

Contact our specialists in pressure ulcer prevention to get more details about the possibilities, uses and effects of the Combi Active Hybrid mattress.



ZiboCare
Denmark

Zibo Combi Active

Comfort and flexibility - with or without a pump



Washable cover

The cover is washable at temperatures up to 95 degrees.

Maximum comfort

Maximum comfort with a viscoelastic top layer on top of the air cells.

Automatic adjustment

Adjusts automatically to produce the optimal therapeutic level for the individual user.

Security border

Specially designed, fixed safety border that minimizes pinching and reduces the danger of falls when getting out of bed.

1 button, 2 treatment options

Without the pump: Self-adjusting pressure distribution.
With pump - alternating pressure levels.

Reduced noise disturbances

Low-noise pump. Easy to use and connect.



Information:

Zibo Combi Active is available in several versions:

80x200x18 cm

83x200x18 cm

85x200x18 cm

88x200x18 cm

Matress

User Weight: 0-250 Kg

Mattress weight: 13 Kg

Composition: 14 air cells.
Top layer: Foam

Cover material: PU
incontinence cover

Cleaning: Soap &
Disinfectants

Pump

Dimensions: 23 x 18 x 8 cm

Pump weight: 1.7 Kg

Settings: 1-2 cycles

Cleaning: Disinfect by wiping with a cloth using common cleaning agents and disinfectants.

Contact us for more information, demonstration and purchase.